

Software Release Notice

Step Proactive v2.7.0

User interface: Step Oncology

Implemented new features and improvements

This section includes the latest news of the new product version.

01. System maintenance mode

The system now supports the manufacturer to put it into a controlled maintenance mode. During this period, app status information is displayed and access is restricted until the system is restored.

Known Issues

This section includes some issues that may occur while using the product and that will be solved in upcoming versions.

Issue 01. It may not be possible to manually enter numerical decimal inputs in Apple brand devices

When a patient tries to report a numerical value including decimals (e.g. fever or weight) with an Apple brand device, decimal separator (e.g. '.' or ',') may not be available on the keyboard.

Workaround:

Use the + or – buttons and the keyboard to obtain the exact decimal number needed.
